



CRISPR Editing of Immortalized Cell Lines with RNPs using Neon Electroporation

Developed by EditCo

Introduction

This protocol describes how to deliver ribonucleoprotein (RNP) complexes that consist of purified Cas9 nuclease duplexed with chemically modified synthetic single guide RNA (sgRNA) to immortalized adherent or suspension cells. An option for knock-in is included. RNP delivery is accomplished using the Thermo Fisher Neon™ Transfection System. A reference for electroporation settings for a wide variety of cell types is included. Chemically modified sgRNAs are designed to resist exonucleases and innate intracellular immune cascades that can lead to cell death. EditCo chemically modified synthetic sgRNAs are of exceptional purity and consistently drive high editing frequencies.

Abbreviations:

CRISPR: clustered regularly interspaced short palindromic repeats

Cas9: CRISPR associated protein 9

sgRNA: single guide RNA

RNP: ribonucleoprotein

PCR: polymerase chain reaction

ICE: inference of CRISPR edits

FACS: fluorescence-activated cell sorting

PBS: phosphate-buffered saline

TE: Tris EDTA

GFP: green fluorescent protein



Materials Required

Material	Ordering Information
Chemically modified sgRNA	EditCo Bio
Cas9 2NLS nuclease (<i>S. pyogenes</i>)	EditCo Bio
Positive control (optional)	Recommended: human <i>RELA</i> sgRNA, <i>CDC42BPB</i> sgRNA (EditCo)
Transfection control (optional)	Recommended: pMAX GFP (Lonza), GFP mRNA (SBI)
TE buffer	Included with EditCo sgRNA
Nuclease-free water	Included with EditCo sgRNA
Neon™ Transfection System	Thermo Fisher Scientific, Catalog #MPK5000
Neon™ Transfection System 10 µl Kit Alternative: Neon™ Transfection System Starter Pack	Thermo Fisher Scientific, Catalog #MPK1025 Thermo Fisher Scientific, Catalog #MPK5000S
Cell counter	Multiple vendors (e.g., Thermo Fisher Scientific)
Normal growth medium	Multiple vendors (cell-type dependent)
1X PBS	Multiple vendors (e.g., Thermo Fisher Scientific)
Tissue culture plates	Multiple vendors (e.g., Thermo Fisher Scientific)
Microcentrifuge tubes	Multiple vendors (e.g., Eppendorf)
TrypLE Express or preferred cell dissociation reagent	Multiple vendors (e.g., Thermo Fisher Scientific)

Important Considerations

Working with RNA and RNPs

- Wearing gloves and using nuclease-free tubes and reagents is recommended in order to avoid RNase contamination.
- Always maintain sterile technique, and use sterile, filter pipette tips.
- All EditCo reagents should be stored according to the manufacturer's recommendations. Synthetic sgRNA should be dissolved in TE buffer and diluted to a working concentration using nuclease-free water. Please consult the [EditCo Gene Knockout Kit Quick Start Guide](#) for best practices related to dissolving and storing synthetic sgRNAs.
- RNP complexes are stable at room temperature for up to 1 hour (may be stored at 4°C for up to one week, or at -20°C for up to 1 month).
Note: that culture medium irrespective of antibiotics. Following transfection, it is not necessary to remove RNP complexes or add/change medium.

Optimized Protocols

- For specific electroporation settings for your cell type, we suggest consulting the [Thermo Fisher Neon™ Transfection System Protocols and Cell Line Data](#).
- Optimization of editing efficiency for a specific cell type will require empirically determining the number of cells required, amount of Cas9 and ratio of sgRNA:Cas9. This guide is meant to provide a starting point for your CRISPR editing experiments.



Suggested Controls

Control	Description	Purpose
Mock	Cells transfected without Cas9 and sgRNA	Wild type sequence for comparison with experimental and other negative controls. Control toxicity from RNP, cell death from electroporation or possible viability issues associated with editing the specific
Negative	Cas9 complexed with a non-targeting sgRNA or no sgRNA	Ensure that there are no false positives due to contamination
Positive	sgRNA with high editing efficiency (e.g., <i>CDC42BPB</i> , <i>RELA</i>)	Ensure all reagents, protocol, and equipment are functioning (effect expected).
Transfection	pMAX GFP (Lonza), GFP mRNA (SBI)	Assess transfection efficiency (without the use of RNPs).

Timeline

Pre-Electroporation		Setup & Electroporation	Post-Electroporation		
Day 1	Day 2	Day 3	Day 4	Day 5	Day 6
Seed Cells Incubate (2 days)		Prepare Destination Plate Assemble RNP Complexes Prepare Cells Transfect cells Incubate (3 days)			Analysis



Protocol

1. Pre-Electroporation

1.1 Seed cells

- a. Subculture cells 2 days before electroporation and seed cells in an appropriately sized vessel so that they are 70-80% confluent on the day of transfection. Each electroporation reaction will require approximately $1 \times 10^5 - 2 \times 10^5$ cells, depending on the cell type.

Note: For cell type specific information, refer to [Thermo Fisher Neon® Transfection System Protocols and Cell Line Data](#).

2. Setup & Lipofection

2.1. Prepare Destination Plate

- a. Pre-warm 1 ml of normal growth medium in each well of a 12-well cell culture plate per reaction. This will serve as the destination plate after electroporation.

2.2 Assemble RNP Complexes (9:1 sgRNA to Cas9 ratio)

EditCo recommends sgRNA:Cas9 ratios between 3:1 and 9:1 for RNP formation. Below is an example using an sgRNA to Cas9 ratio of 9:1 for a single reaction (scale up appropriately).

- a. Prepare sgRNA stock at 30 μM and Cas9 nuclease stock at 20 μM , and store at -80°C until use.
- b. In appropriate plates/ tubes, assemble RNP complexes in the order shown below.
- c. Incubate RNPs for 10 minutes at room temperature.

Note: The sgRNA:Cas9 ratios may need to be determined empirically to achieve optimal editing efficiency.

RNP Components, Molarity, & Volume		
Component	Molarity	Volume (per reaction)
sgRNA	30 μM (pmol/ μl)	3 μl (90 pmol)
RNP Components, Molarity, & Volume	20 μM (pmol/ μl)	0.5 μl (10 pmol)
Cas9	3 μM (pmol/ μl)	3.5 μl
Resuspension buffer	-	7 μl
Total volume	-	28.3 μl

Knock-in Option: to knock in small inserts (<50 bp), an ssDNA HDR Template can be added. The recommended length of each homology arm is at least 50 bp. Add 1 μl 60 μM ssDNA HDR Template per reaction to each well. Optimization may be required.



2.3 Prepare Cells

Note: For suspension cells: spin down cells before each aspiration of culture medium and washes (step a below). Skip steps b and c below.

- Aspirate cell culture medium and wash cells 1-2 times with 1X PBS.
 - Add TrypLE Express and incubate the cells for ~5 minutes, or until they detach from the plate completely.
- Note:** Do not shake or hit the flask to dislodge cells, as this may lead to clumping and inaccuracies in cell counting.
- Neutralize the dissociation reaction with 2X volume of normal growth medium.
 - Count cells to determine the cell density (using cell counter).
 - Transfer 1-2 x 10⁶ cells to a sterile microfuge tube. One tube will contain enough cells for ~10 transfections.
 - Centrifuge cells for 5 minutes at 500 x g. Aspirate medium.
 - Wash the cells once with 1X PBS and repeat the centrifugation step. Aspirate PBS.
 - Resuspend the cell pellet in 50 µl of resuspension buffer R (provided with Neon™ Transfection System 10 µl Kit).

Note: Avoid storing the cell suspension for more than 15 minutes at room temperature, as this reduces cell viability and transfection efficiency.

- Add 5 µl of cell suspension to each RNP solution (7 µl) to make 12 µl of cell-RNP solution per reaction.

2.4 Transfect Cells

- Aspirate 10 µl of cell-RNP solution to a 10 µl Neon tip.
- Electroporate using cell type optimized conditions.

Note: Refer to [Thermo Fisher Neon® Transfection System Protocols and Cell Line Data](#).

- Immediately transfer cells to a pre-warmed 12-well plate (prepared in step 2.1)
- Incubate the cells for 2-3 days in a humidified 37°C/5% CO₂ incubator.

3. Post-Lipofection

2.5 Analysis

- Extract DNA from cells.
- Conduct analyses to determine editing efficiency: PCR, Sanger sequencing, and [ICE analysis](#). Next-Gen Sequencing, FACS, or functional tests may be conducted as alternatives.

Option: If storing cells for future use is desired, split cells into two groups (one for analysis and one for cell culture).

Additional Information

For an up-to-date list of all protocols and other resources, please visit this [link](#).

For technical assistance, contact our Scientific Support Team at technicalsupport@editco.bio.

For common FAQs, please visit this [link](#).